USER MANUAL




1) General Remarks

- It is essential that you and any other operator of this product read and understand the content of this manual before installing and using this product !
- In order to minimize operation errors, this manual must be access able to the staff at all times.
- The DURANBOARD is only to be used as instructed.
- Please refer particularly to the user manual of the transport-stretcher.
- · Descriptions do not necessarily depict contents and are not true to scale.
- We take no liability for damages caused by operating errors or incorrect assembly.
- Please pay close attention to the country-related, applicable safety regulations for patient transfer.
- Subject to technical changes.



2) Technical Data

- 1... Pneumatic suspensions with automatic weight adjustment from 40-250 kg and hydraulic shock absorber. Stroke distance 95 mm.
- 2... tailgate with cylindrical roller
- 3... Loading flap
- 4... Switch box, contains: a) compressor 12V DC, 16A with check valve and integrated thermal overload protector.
 - b) relay 12V DC
 - c) magnetic valves 12V DC
 - d) terminal block
 - e) noise insulation
- 5... Air pressure tank (4 liters) with pressure switch
- 6... Stretcher mounting
- 7... Base frame
- 8... Gas spring





2) Technical Data

- General-purpose system for assembly of various transport-stretchers (e.g. Stryker, Stollenwerk, Ferno aso.) combined with the stretcher's fastening system.
- System height: 140 mm (lowered) 250 mm operating height
- The total height varies depending on vehicle specifications
- Total length including closed loading ramp 2130mm
- Total weight approx. 80kg
- Maximum loading capacity approx. 400kg (incl. stretcher)
- Ignition AND main switch ON: Device ready
 - Ignition OR main switch OFF: Device lowered (eg. for reanimation)
- Electric connection: brown = ground (wire no. 31) red = positive (wire no. 15)
- Current demand max. 20A at 12 Volt DC
- Control circuit fused 5A

3) Disinfection, Cleaning and Maintenance

- Disinfection and cleaning of the DURANBOARD are to be performed according to the regulations of the respective country.
- The DURANBOARD does not require further maintenance. In case of faults or queries, please consult our hotline.
- For complete please refer to repair instructions to download under: <u>www.duranboard.com</u>



4) Starting Up

- After successful assembly and electric connection, start ignition of vehicle and turn on the main switch (9).
- The compressor will fill the air supply tank (5) to approx. 8,2 bar within about 2 minutes and the DURANBAORD is ready.
- The DURANBOARD works automatically when the air pressure has lowered to approx. 7,5 bar, the compressor starts automatically for approx. 15 seconds again. This allows constant operating pressure.
- By opening the loading ramp (2) and/or turning off the ignition, the DURANBOARD is automatically lowered to allow an easy and energy-efficient loading and unloading.
- For correct loading and unloading please refer to the instructions of the (transport-) stretcher.
- Pay attention to the correct locking of the (transport-) stretcher on the DURANBOARD.
- For reanimation, the main switch (9) is to be turned off, the DURANBOARD will lower immediately.

ATTENTION!

As a precaution, the stretcher should always be held when opening the loading ramp !!




Adress: A-6130, Schwaz, Münchner Strasse 30

Hotline: +43 676 928 71 07 Email: info@duranboard.com



5) Circuit Diagram




6) Connection Layout

target	colour	cross-section		target	colour	cross-section
	yellow blue	1,0	0 0 0 0	magnetic valve magnetic valve relay (85)		1,0 1,0 1,0
			0 0	pressure switch pressure switch	grey grey	1,0 1,0
	yellow black	1,0 1,0	• • • • • •	main switch magnetic valve magnetic valve		1,0 1,0 1,0
jumper	black	1,0	0 0	compressor ground	black brown	2,5 2,5
	orange red	1,0 2,5	• • • • • • • • • • • • • • • • • • •	fuse 5A	red	1,0
relay (87)	red	2,5	0 0	fuse 5A	red	1,0
relay (30)	red	2,5	0 0	positive 25A	red	2,5

Hotline: +43 676 928 71 07 Email: info@duranboard.com



7) Pneumatic Connection



tube colours

Red	 working circuit 	
Blue	= inter-valve-circui	t
Black	 supply circuit 	

Parts

- A... compressor
- B... air supply tank C... pressure switch
- D... level control valve
- E... air cushion
- F... magnetic valve
- G... check valve
- P ... pressure supply
- 1... magnetic valve > supply
- 2... magnetic valve > working press.
- 3... air release



8) Spare Parts



30205 compressor set

30063 pressure switch, pre-finished

30054 magnetic valve, pre-finished



30053 lock valve, pre-finished



30291 control valve



30054 connection valve



30263 main switch 30251 fuse holder 30252 blade-type fuse 25A 30253 blade-type fuse 5A



Adress: A-6130, Schwaz, Münchner Strasse 30

Hotline: +43 676 928 71 07 Email: info@duranboard.com DB

8) Spare Parts



30700 swiveling system 30081 gaiter 30311 air cushion



30053 vibration damper 30110

gas compression spring



brass connector



30342 cable gland 30343 fence jack 50150 extension carrier

Adress: A-6130, Schwaz, Münchner Strasse 30

Hotline: +43 676 928 71 07 Email: info@duranboard.com



8) Spare Parts



30220 noise insulation box



rubber buffer

30500

connector 2 wires



30170 terminal block



30150 wire protective tube



30350 plug 30 x 30 30361 plug 60 x 40



30240 pressure hose red 30241 pressure hose blue 30242 pressure hose black



30243 thermal protective tube

Adress: A-6130, Schwaz, Münchner Strasse 30

Hotline: +43 676 928 71 07 info@duranboard.com Email:

